

**Clean Version of Amended Paragraph beginning on Page 32, Line 32**

---

C/ The selected material 170 can also comprise demineralized bone matrix suspended in glycerol (e.g., Grafton<sup>TM</sup> allograft material available from Osteotech), or SRS<sup>TM</sup> calcium phosphate cement available from Norian. These viscous materials, like the bone cement previously described, can be loaded into the syringe 104 and injected into the cavity using the nozzle 106, which is inserted through the cannula instrument 30 into the cavity. The tamping instrument 108 is used to displace residual material from the cannula instrument 30 into the cavity, as before described.

---